

PCR

Introduction

Polymerase chain reaction (**PCR**) is a laboratory technique used to make multiple copies of a segment of DNA. **PCR** is very precise and can be used to amplify, or copy, a specific DNA target from a mixture of DNA molecules.

Materials

- › Phusion DNA Polymerase, 2 U/μL
- › 5X Phusion HF Buffer
- › 5X Phusion GC Buffer
- › 50 mM MgCl₂ solution
- › DMSO
- ›

Procedure

1. CR reactions should be set up on ice.
2. Prepare a master mix for the appropriate number of samples to be amplified.
3. The Phusion™ High-Fidelity DNA Polymerase should be pipetted carefully and gently as the high glycerol content (50%) in the storage buffer may otherwise lead to pipetting errors.
4. Due to the nature of the Phusion™ High-Fidelity DNA Polymerase, the optimal reaction conditions may differ from PCR protocols for standard DNA polymerases.
5. Due to the high salt concentration in the reaction buffer, the Phusion™ High-Fidelity DNA Polymerase tends to work better at elevated denaturation and annealing temperatures.

	A	B	C	D
1	Component	20 μL rxn	50 μL rxn	Final conc.
2	H ₂ O	add to 20 μL	10 μL	1X
3	5X Phusion™ HF Buffer	4 μL	1 μL	200 μM each
4	Forward primer	X μL	X μL	0.5 μM
5	Reverse primer	X μL	X μL	0.5 μM
6	Template DNA	X μL	X μL	
7	(DMSO, optional)	(0.6 μL)	(1.5 μL)	(3%)
8	Phusion™ High-Fidelity DNA Polymerase	0.2 μL	0.5 μL	0.02 U/μL
9	10 mM dNTPs	0.4 μL	1 μL	200 μM each

- 6.

	A	B	C	D	E	F
1	Cycle step	2-step protocol		3-step protocol		Cycles
2		Temp.	Time	Temp.	Time	
3	Initial Denaturation	98°C	30 s	98°C	30 s	1
4	Denaturation	98°C	5–10 s	98°C	5–10 s	25-35
5	Annealing	-	-	X°C	10–30 s	
6	Extension	72°C	15–30 s/kb	72°C	15–30 s/kb	
7	Final extension	72°C	5–10 min	72°C	5–10 min	1
8	Hold	4°C	Hold	4°C	Hold	Hold

7.