

# Site Directed Mutagenesis

---

## Introduction

The Q5 Site-Directed Mutagenesis Kit enables rapid, site-specific mutagenesis of double-stranded plasmid DNA in less than 2 hours (Figure 1). The kit utilizes the robust Q5 Hot Start High-Fidelity DNA Polymerase along with custom mutagenic primers to create insertions, deletions and substitutions in a wide variety of plasmids. After PCR, the amplified material is added directly to a unique Kinase-Ligase-DpnI (KLD) enzyme mix for rapid (5 minutes), room temperature circularization and template removal. Transformation into high-efficiency NEB 5-alpha Competent E. coli, provided with the kit, ensures robust results with plasmids up to at least 14 kb in length.

## Materials

- › Q5 Hot Start High-Fidelity 2X Master Mix (2X)
- › KLD Enzyme Mix (10X)
- › KLD Reaction Buffer (2X)
- › Control SDM Primer Mix (10  $\mu$ M each)
- › Control SDM Plasmid (5  $\mu$ g/ml)
- › NEB 5-alpha Competent E. coli (High Efficiency)
- › pUC19 Transformation Control Plasmid (50 pg/ $\mu$ l)
- › SOC Outgrowth Medium
- ›

## Procedure

### Step I: Exponential Amplification

1. Assemble the following reagents in a thin-walled PCR tube.

	<b>A</b>	<b>B</b>	<b>C</b>
1		25 $\mu$ l RXN	FINAL CONC
2	Q5 Hot Start High-Fidelity 2X Master Mix	12.5 $\mu$ l	1X
3	10 $\mu$ M Forward Primer	1.25 $\mu$ l	0.5 $\mu$ M
4	10 $\mu$ M Reverse Primer	1.25 $\mu$ l	0.5 $\mu$ M
5	Template DNA (1–25 ng/ $\mu$ l)	1 $\mu$ l	1-25 ng
6	Nuclease-free water	9.0 $\mu$ l	

2. Mix reagents completely, then transfer to a thermocycler.

3. Perform the following cycling conditions:

	<b>A</b>	<b>B</b>	<b>C</b>
1	STEP	TEMP	TIME
2	Initial Denaturation	98°C	30 seconds
3	25 Cycles	98°C	10 seconds
4		50–72°C*	10–30 seconds
5		72°C	20–30 seconds/kb
6	Final Extension	72°C	2 minutes
7	Hold	4–10°C	

## Step II: Kinase, Ligase & DpnI (KLD) Treatment

4. Assemble the following reagents:

	<b>A</b>	<b>B</b>	<b>C</b>
1		Volume	Final Conc
2	PCR Product	1 $\mu$ l	
3	2X KLD Reaction Buffer	5 $\mu$ l	1X
4	10X KLD Enzyme Mix	1 $\mu$ l	1X
5	Nuclease-free Water	3 $\mu$ l	

5. Mix well by pipetting up and down and incubate at room temperature for 5 minutes.

## Step III: Transformation

6. Thaw a tube of NEB 5-alpha Competent E. coli cells on ice.

7. Add 5  $\mu$ l of the KLD mix from Step II to the tube of thawed cells. Carefully flick the tube 4-5 times to mix. Do not vortex.

8. Place the mixture on ice for 30 minutes.

9. Heat shock at 42°C for 30 seconds.

10. Place on ice for 5 minutes.

11. Pipette 950  $\mu$ l of room temperature SOC into the mixture.

12. Incubate at 37°C for 60 minutes with shaking (250 rpm).

13. Mix the cells thoroughly by flicking the tube and inverting, then spread 50-100  $\mu$ l onto a selection plate and incubate overnight at 37°C. It may be necessary (particularly for simple substitution and deletion experiments) to make a 10- to 40-fold dilution of the transformation mix in SOC prior to plating, to avoid a lawn of colonies.