

Anti-Invertase Assay

Project: Sequencing

Authors: Suhas Badadal

Created at: 2020-07-09T14:08:01.665323+00:00

THURSDAY, 9/7/2020

Aim:

To verify the action of anti-invertase on invertase, and determine the concentration needed to inhibit invertase activity.

Materials Required:

- 96 well plate
- 100 ml conical flask
- LB media
- Cloning kit
- Glucose stock
- Sucrose stock
- Pipettes
- 10 μ l microtips
- 100 μ l microtips

Procedure:

- Clone anti-invertase-histag gene into E.coli K-12 using pET28-a backbone as directed by the [Cloning](#) protocol, and plates are prepared (Kanamycin)
- Prepare culture using LB media as directed by the [LB Medium Preparation](#) protocol
- Anti-invertase protein is extracted as directed by the [His-tag Antibody protein purification] protocol and stored at -20°C . (1 L culture is expected to give around 2-3 mg, 4 mg at best, while a 10 L culture can give 25 mg.)
- A commercial **invertase assay kit** employing yeast invertase (Cat.MAK118, Sigma Aldrich, USA) is used as per the manufacturer's protocol to examine the potency of recombinant ShINH1.
<https://www.sigmaaldrich.com/catalog/product/sigma/mak118?lang=en®ion=IN>
- 40 μ l of reaction volume containing different concentrations of ShINH1, from 0.1 μM to 0.5 μM , in increments of 0.1 μM are pre-incubated in a 96-well plate with commercial acid invertase at 37°C for 30 min.
- Glucose standards (40 μ l volume, 0–100 μM glucose) are added to separate wells of the plate. The same volume of reaction buffer was used as the assay blank in separate wells.
- Substrate was added to each well 5 μ l of 20 mM sucrose followed by incubation for 20 min at room temperature.
- After incubation, the reaction mixture containing 95 μ l of reaction buffer, 1 μ l of enzyme mix and 1 μ l of dye reagent (all supplied with the kit) was prepared and 90 μ l of reaction mixture was added to each of the blank, sample, and standard wells followed by incubation for 20 min at room temperature in darkness.
- The amount of glucose liberated was calculated from the glucose standard curve.
- The specific activity of enzyme was calculated and expressed as μmoles of glucose formed per milligram of protein per minute.

Procedure							
	A	B	C	D	E	F	G
1	Sample	[Anti-Invertase Protein] (in μM)	Reaction volume of Anti-Inv (in μl)	[Sucrose] (in mM)	Volume of [Sucrose] stock to be added (in μl)	Volume of reaction mix to be added (in μl)	Total Volume (in ml)
2	Blank	0.00	40	20.00	5.00	90.00	0.135
3	1	0.10	40	20.00	5.00	90.00	0.135
4	2	0.20	40	20.00	5.00	90.00	0.135
5	3	0.30	40	20.00	5.00	90.00	0.135
6	4	0.40	40	20.00	5.00	90.00	0.135
7	5	0.50	40	20.00	5.00	90.00	0.135
8	Total	1.50	240	120.00	30.00	540.00	0.81